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EN User manual

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EN: BrandsLAB sp. z o.o reserves the right to make changes and improvements to its product without prior notice. Visit our website for the latest version of the manual.

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www.brandslab.pl/productsupport

For assistance or warranty claims, please contact us using the form on our website

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MADE IN CHINA

EN

1. Product Description and Intended Use

Flexible Hot & Cold Gel Pack the gel pack conforms to the shape of the body and is intended for home use.

The product is designed to provide a cooling or warming sensation to selected areas of the body.

2. Safety Information and Precautions

Please read these instructions before first use and keep them for future reference.

- For external use only.
- Do not puncture or cut.
- Do not use if damaged or leaking.
- Do not use on children without adult supervision.
- Do not use while sleeping.
- Check the product temperature before use.
- Do not overheat or exceed a total heating time of 20 seconds.
- Do not freeze. Cool only in a refrigerator.
- Discontinue use immediately if discomfort or irritation occurs.
- Keep out of reach of children.

3. Directions for Use

Cooling

1. Place the gel pack in the resealable bag provided to keep it clean and hygienic.
2. Cool in a refrigerator for 1–2 hours.
3. Remove from the bag and apply to the desired area of the body for up to 20 minutes.

Note: Do not freeze. Cooling in a refrigerator ensures proper product performance.

Heating

1. Heat in a microwave oven at medium power for 15–20 seconds.
2. If additional warmth is required, heat in additional 10-second intervals.

3. Do not exceed a total heating time of 20 seconds.
4. Check the product temperature before use.
5. Apply for 15–20 minutes.
6. Allow the product to return to room temperature before reheating.

4. Cleaning and Storage

Composition:

Outer layer: 80% Nylon, 20% Spandex

Filling: Hydrogel (55% Water, 45% Polyacrylamide)

- Do not wash.
- Do not bleach.
- Do not dry clean.
- Do not tumble dry.
- Do not iron.

Clean with a damp cloth and mild detergent.

Store in a dry, well-ventilated place away from heat sources and direct sunlight.

5. Disposal

The product is made from multiple materials (textiles and hydrogel)

Do not dispose of the contents into drains or waterways.

Dispose of the used or damaged product at a municipal waste collection facility or in accordance with local waste segregation requirements applicable at the place of use.

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